



JANNY MT



Flexible Installations of Controlled Atmosphere Units

Janny MT 610L Module



Janny MT 610L Module

“Janny MT preservation”

- “**Janny MT Module**” consists in a Tiempo Cap and a plastic pallet-box of **610 liters**. The unit forms an individual module with controlled atmosphere. This can be used for long storage of fruits, flowers, vegetables or mushrooms.
- Thanks to the passive diffusion of gases through the membrane, **Janny MT Module** provides an optimal atmosphere meeting the storage requirements of the chosen produce.

Characteristics

- Patented closing by clipping
- Cap equipped with membranes
- Easy to use
- Full rigid box: 610 liters (strengthened industrial box).
- Payload: 600 kg/1322.77 lb
- Stacking: up to 2000 kg/4409.25lb on cap.
- Washable
- Can be used for up to 15 years
- Food quality certified
- ABS cap
- HDPE Box
- No metal component
- Water tightness ensured by silicon seal

Use

- **Janny MT Module** can be used to store fruits, vegetables, mushrooms and flowers.

Product data	Unit	Value
External dimensions	cm	101 x 121 x 79
	inch	39.76 x 47.64 x 31.1
Internal dimensions	cm	92 x 112 x 56,5
	inch	36.22 x 44.1 x 22.24
Capacity	liter	610
Hold on the ground	feet	4

Application characteristics

Tiempo Cap must be installed after cooling the produce. So as to maintain the characteristics of the atmosphere created, it is not recommended to open the box during storage time. It is essential to monitor the evolution of the atmosphere during storage time by measuring the O2 and CO2 levels. This can be done with the Tiempo Test device

Precaution for use

When setting up, make sure the proper fastening the lid on the container to seal. If JANNY MT module is not used in recommended conditions, product performance may be modified. The sealing cover is essential for the passive diffusion of gases. Particular attention should be paid to the joint membranes and during the non-use of the covers. The JANNY MT module is used for storage of live products (breathing).

Sté Janny MT SAS
La Condemine 71260 PERONNE
FRANCE

Tel: +33.3.85.23.96.20

Fax: +33.3.85.36.96.58

www.jannymt.com

jannymt@jannymt.com

@Janny.MT.SAS